

How to Wear your Patti Mask (S/M/L)

How you wear your mask affects its effectiveness!

The Patti Mask is a **Halyard Mask** available in small/medium or large size. It is adapted from the design of UF Health Prototype 2 **modelled after a N95 mask**. It is **reusable** and **as effective as a surgical mask for filtering particles** and **as effective as a N95 mask for creating a tight seal** around the nose and mouth. The cup shape is to **maximize breathability** by maximizing the surface area for easy breathing.





1. Shape the metal nose piece to create a tight seal around your nose

Flame Resistance











2. Tie the top straps to secure them around your head just above the ears

3. Tie the bottom straps to secure them around your head above the first tie



4. Pull the bottom of the mask slightly over your chin

NIOSH PFE 60-80%

5. Check your seal: take a deep breath in and out through your mouth. The mask should slightly collapse inward when you breathe in and re-expand when you breathe out

You can request your free mask by filling out a form at freemaskprojectvancouver.com/donate-now

DISCLAIMER: THE MASK INFORMATION IS PROVIDED "AS-IS, WHERE-IS," WITHOUT REPRESENTATIONS, CONDITIONS OR WARRANTIES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE MASK INFORMATION WILL NOT INFRINGE ANY PATENT, COPYRIGHT, TRADEMARK OR OTHER PROPRIETARY RIGHTS. THE RECIPIENT IS SOLELY RESPONSIBLE FOR DETERMINING THE APPROPRIATENESS OF USING, REPLICATING OR REDISTRIBUTING THE MASK INFORMATION AND DESIGN. IN THIS REGARD, THE RECIPIENT ASSUMES ALL LIABILITY FOR DAMAGES, OF WHATEVER NATURE AND DESCRIPTION, WHETHER IN CONTRACT OR IN TORT, WHICH MAY ARISE FROM THE USE OF THE MASK INFORMATION AND DESIGN. THE FREE MASK PROJECT OF VANCOUVER, INCLUDING ITS EMPLOYEES AND AGENTS, WILL NOT BE LIABLETO THE RECIPIENT OR TO ANY THIRD PARTY FOR ANY LOSS, CLAIM OR DEMAND MAD DESIGN BY THE RECIPIENT, OR ANY LOSS,

The Free Mask Project of Vancouver

freemaskprojectvancouver.com